

Sung Kwang Meditech

For surgical patient & Rehabilitation Radiology Products

Operating Room


Snow pac (Bean Bag) Several device Position pad Gel product Instrument Case Silicon product



Rehabilitation Gel wheelchair cushion pad Gel wheelchair back cushion pad

Silver Gel embo cushion mattress





Radiation Protection For: Cardiac Cath Lab Interventional Radiology Pain Management









T) 82-31-855-7548 F) 82-31-855-7549 E) skmeditech@skmeditech.com A) 1014-20, Bukhansan-ro, Jangheung-myeon, Yangju-si, Gyeonggi-do, Korea



www.skmeditech.com


SNOW PAC (Bean Bag)

FOR SURGERY POSITIONING PATIENTS

Snow Pac are soft pads (Available in various sizes and shapes) that are filed with thousands of tiny plastic beads. Snow Pac provides more stable and comfortable position to patients and saves valuable time for the entire O.R team on every case. It prevents pressure points (With Gel). It aids faster, easier set-up and repositioning during surgery.

It aids to become comfortable position. It aids better surgery position presentation.

PRINCIPLE OF SYSTEM





EFFECT

Stable

Because Snow-Pacs are molded well, firmly and securely to shape of the patient's body than sandbags, towel rolls and other devices that are merely propped against the patient.

Comfortable

The apparent wrinkles are inverted-drawn inward by suction and will not irritate the skin.

Better Presentation

With Snow-Pacs the patiant can be positioned exactly as the surgeon desires for excellent presentation of the operative site. Snow-Pacs conform to any shape, so the patient's position does not have to be compromised.

More Versatile

With Snow Pacs patient can often be placed in positions that could not be easily achieved using other positioning devices.

No Pressure Points (With Gels)

A Snow-Pac supports all parts of the patient's body uniformly when used with gel like a floating mattress.

Saves Time

Patients can be positioned much more quickly and easily using Snow Pacs, and if necessary, be repositioned during a case.

Before

Fix or put the Snow Pac to the patient's body according to the position you want to form (e.g., Supine Position) and connect the suction line to the entrance of the Snow Pac.

Depending on the intensity of the Snow Pac you want to make, use a suction line to inhale the air in the Snow Pac properly. (If you inhale for about 10 to 15 seconds, vaccum is created inside the Snow Pac and it becomes hardens)

After

The inside of the Snow Pac is vacuumed and fixed to the desired shape.

When the suction line is disconnected, the Snow-Pac retains its firm, molded shape and holds the patient securely in position.

After surgery, if user open the entrance of the Snow Pac and press it slightly, air re-enters into the Snow Pac and it becomes soft again and re-useable.

OTHER FEATURES

Any Vacuum Source

Use wall suction, aspirator pump or any other convenient O.R. vacuum source.

Reusable

As the exterior material is soft and sturdy, Snow Pac can be used for many years as long as they are thoroughly controlled for contamination.

Radiolucent

Excellent radio transparency with only minimal and clearly identifiable artifacts.

Clean, Disinfect, Sterilize

Can be cleaned with regular soaps and liquid disinfectants or sterilized by gas.

Saves Money


Snow Pacs save the cost of sandbags, towel rolls and other positioning devices, and they save time in every case

Repair Service

Factory inspection and repair service is available at nominal cost after warranty period. Full One Year Warranty.


INDICATION

KIDNEY POSITION, LATERAL POSITION



THYROID POSITION





HEAD STABILIZATION POSITION





SUPINE POSITION



VARIOUS FULL BODY POSITION



Sung Kwang Meditech Co.,Ltd	www.skmeditech.com	Contact to E-mail: skmeditech@skmeditech.com	Sung Kwang Meditech Co.,Ltd	www.skmeditech.
		Tel: 82-31-855-7548 Fax: 82-31-855-7549		

Sung Kwang Meditech $C \in$



A-1-2 type A-1-3 type

Easily up in nephrectomy (A type), and also for handicapped person with cerebral palsy.(A-1 type).



A type	73(L) * 100(W)cm	
A-1 type	73(L) * 100(W)cm	
A-1-1 type	90(L) * 100(W)cm	
A-1-2 type	73(L) * 100(W)cm(U)	
A-1-3 type	90(L) * 100(W)cm(U)	

Provides hyperextension, prevents postoperative discomfort.

B type	72(L) * 88(W)cm	
B-1 type	118(L) * 88(W)cm	

Stabilizes head in eye, facial or dental surgery.

30(L) * 40(W)cm C type



Supports and eliminates back pains during long procedures. Prevents foot drop splint in

Emergency.

30(L) * 55(W)cm D type

E type

E type 165(L) * 72(W)cm

POSITION PADS



PAD FOR PRONE POSITION




Chest pad (t	hick, thin)
SK-5-1	18(L) * 52(W) * 13(H)cm
Head pad	
SK-6	32(L) * 25(W) * 10(H)cm
SK-6-1	24(L) * 26.5(W) * 18(H)cm
Knee pad	
SK-7	30(L) * 20(W) * 6(H)cm
Leg pad	
SK-3	40(L) * 17(W) * 17(H)cm

Contact to E-mail: skmeditech@skmeditech.com

Tel: 82-31-855-7548 Fax: 82-31-855-7549

SK-6-1



SK-3

Sung Kwang Meditech Co., Ltd

www.skmeditech.com

Good positioner Anti-bedsore effect Very soft, Good elasticity Good endure.

Poly urethane vinyl(Outer) Water resistant.

Special sponge(Inner)

SUPINE	
Head, foot pad	
SK-1	22(L) * 16(W) * 11 (5.5)(H)cm
Head pad	
SK-1-1	22.8(0.D) * 11.5(I.D) * 5(H)cm
Sh oulder pad	
SK-2	23(L) * 14(W) * 7(H)cm



Instrument Sterilization Cases are designed to protect fine surgical instruments and various scopes during sterilization and storage. Our cases are made of elastic special materials for flexible handling and it is transparent so that the inside can be seen during autoclaving. There are fine circular holes in the upper(Lid) and lower(Base) of the case to help circulate steam during sterilization. Silicon Embo Mat or Bar is placed on the bottom and prevents instruments from moving freely.

- This sterilizing instrument tray is designed to withstand repeated exposure to 134 degrees (Steam / EO / Plasma steralizable).
- There are a variety of designs and sizes suitable for device applications such as micro instrument, Laparoscope, Telescope, and Endoscope.
- The size table below is based on the outer diameter, and the inner diameter by -1 to 2cm.

SK-19(2F)	$26(L) \times 17(W) \times 4(H) \text{ cm}$	
SK-19-a(1F)	26(L) x 17(W) x 2.5(H) cm	
	₩ 100 000000000000000000000000000000000	

Micro Instrument Case(1F)		
SK-19-S	16(L) × 6(W) × 3(H) cm	
SK-19-1	25(L) x 5(W) x 2.5(H) cm	
SK-19-3	30(L) x 10(W) x 2.5(H) cm	




SK-21-1

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com


Sung Kwang Meditech $C \in$



CASES

INSTRUMENT STERILIZATION TRAYS USE CASES

The instrument sterilization case photos below are examples of the use of cases that were actually applied in storage or sterilization processes by various departments such as surgery and internal medicine.



SK-19-5(2F) 39(L) x 27(W) x 4(H) cm



37(L) x 15(W) x 7(H) cm



SK-19-a(1F) 26(L) x 17(W) x 2.5(H) cm

Strengths of Sungkwang Meditech trays:

- Cases are transparent. (Users can see instruments during autoclaving)
- Cases are light but strong. (We do not use stainless, steel materials. But, use special strong materials)
- There are fine circular holes in the upper(Lid) and lower(Base). (It helps circulate steam during sterilization)
- Silicon Embo Mat or Bar is available.
- (It is placed on the bottom and prevents instruments from moving freely). • It is possible to assemble freely based on the shape of the instrument.





SK-20-3 37(L) x 15(W) x 7(H) cm

SK-21-1 64(L) x 25(W) x 15.5(H) cm









SK-20-3 37(L) x 15(W) x 7(H) cm



Contact to E-mail: skmeditech@skmeditech.com Sung Kwang Meditech Co.,Ltd Sung Kwang Meditech Co.,Ltd www.skmeditech.com www.skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549 6



SK-19-5-b(1F) 39(L) x 27(W) x 12(H) cm



SK-21-1 64(L) x 25(W) x 15.5(H) cm



SK-19-3 30(L) x 10(W) x 2.5(H) cm


Sung Kwang Meditech $C \in$



SK-20-5

SK-20-3 37(L) x 15(W) x 7(H) cm

SK-20-4 $60(L) \times 20(W) \times 6(H) \text{ cm}$



SK-20-4 $60(L) \times 20(W) \times 6(H)$ cm


CASES

NEEDLE COUNTER & SHARPS COLLECTOR

Inner : magnet For surgery in OR It aids to count needle thoroughly.



It finds needles by rolling on the floor



SK-18 (S)	14(L) * 4(W) * 1.5(H)cm	
SK-18-1(L)	17(L) *8(W) *2(H)cm	
SK-13-2	100(L) * 30(W)	

MICRO INSTRUMENT ONE FLOOR CASE





MICRO INSTRUMENT ONE FLOOR CASE

SK-18





SK-18-1

SK-19-a (1F)

(External) SK-19-a(1F) 26(L) * 17(W) * 2.5(H)cm

26(L) * 17(W) * 4(H)cm

(External)

MICRO INSTRUMENT TWO FLOOR CASE





Contact to E-mail: skmeditech@skmeditech.com Sung Kwang Meditech Co.,Ltd www.skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549 Sung Kwang Meditech Co.,Ltd

www.skmeditech.co

8

Sung Kwang Meditech $\zeta \in$

For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable




SK-19-3


(External)

SK-19-S	16(L) * 6(W) * 3(H)cm
SK-19-1	25(L) * 5(W) * 2.5(H)cm
SK-19-3	30(L) * 10(W) * 2.5(H)cm

com	Contact to E-mail: ski	meditech@skmeditech.cc	m
	Tel: 82-31-855-7548	Fax: 82-31-855-7549	9

MICRO INSTRUMENT ONE FLOOR CASE







MICRO INSTRUMENT TWO FLOOR CASE

For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

(External


Steam / EO / Plasma Sterilizable

SK-19-5-a(1F) 38(L) * 26(W) * 2.5(H)cm SK-19-5-b(1F) 39(L) * 27(W) * 12(H)cm

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

LAPAROSCOPE CASE



TELESCOPE CASE





Steam / EO / Plasma Sterilizable

(External) 39(L) * 27(W) * 4(H)cm SK-19-5(2F)

Tel: 82-31-855-7548 Fax: 82-31-855-7549

Contact to E-mail: skmeditech@skmeditech.com

Sung Kwang Meditech Co.,Ltd

SK-19-5 (2F)

www.skmeditech.com

Sung Kwang Meditech Co.,Ltd


www.skmeditech.com

10

Sung Kwang Meditech $\zeta \in$

SK-20-1

For Laparoscope. Silicon bar is placed to prevent



For Telescope.

For Arthroscope Silicon bar is placed to prevent laparoscope to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable




(External)

SK-20-2

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549 11

28(L) * 7(W) * 5(H)cm

ENDOSCOPE INSTRUMENT CASE




LONG SCOPE INSTRUMENT CASE





LAPAROSCOPIC INSTRUMENT CASE

Sung Kwang Meditech Co.,Ltd



www.skmeditech.com

Tel: 82-31-855-7548 Fax: 82-31-855-7549

12

Sung Kwang Meditech $C \in$

Sterilization cases used during laparoscopic long scope and robotic instrument surgery.

Silicon bar or mat is placed to prevent movement of instruments used in obstetrics, urology, otorhinolaryngology, surgery and robotic surgery.

Two or more laparoscopic instruments can be inserted.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable


(External)

SK-20-5 72(L) * 15(W) * 10(H)cm

Silicon bar is placed to

	(External)
SK-21-1	64(L) * 25(W) * 15.5(H)cm
SK-21-1-A	64(L) * 25(W) * 12.5(H)cm

OPHTHALMIC INSTRUMENT CASE



For microsurgical and other delicate instruments, for maximum protection and flexible handling, made of resilient polymer.

Protruded silicon embomat is placed on the bottom prevents instrument to move freely.

Cases are designed to withstand repeated exposure to 134°C steam temperature.

Steam / EO / Plasma Sterilizable

(External)

SK-19-S 16(L) * 6(W) * 3(H) cm SK-19-a (1F) 26(L) * 17(W) * 2.5(H)cm SK-19-5-a (1F) 38(L) * 26(W) * 2.5(H)cm

Material : stainless steel with marked L/N number, Use : for specimen storage

order made available

BIOPSY CASE

Delphain	Pyramidal	L2
		L3
R. Thyroid	L. Thyroid	L. Jungular
		L4
R.CN	L.CN	L5
		L. Para. T
	R. Thyroid	R. Thyroid L. Thyroid

SK-21-a-1 (for Thyroid)





2 3 4 5 6 7 4sa

13 14a 14v 16 16a2 16b1 17

8 9 10 11 11p 11d 12 14v



4sb 4d

SK-21-a-3 (for stomach)

12a 12b 12c 12h 12p

SK-21-a-4 (for biliary and Pancreas)

(External) SK-21-a-1 for Thyroid 16.5(W) * 19(L) * 3.2(H)

SK-21-a-2 for breast 14.5(W) * 15(L) * 3.2(H) SK-21-a-3 for stomach 31.5(W) * 8.5(L) * 3.2(H) SK-21-a-4 for biliary and Pancreas 30(W) * 8.5(L) * 3.2(H)

Operating Room

GEL PRODUCTS

Capillary closure is at the root of ulcer formation.

The blockage or damage of capillaries prevents the natural maintenance of healthy skin. In operation, the use of proper positioning devices prevent excessive pressure build-up and capillary closure. Evenly distribute weight and provide constant support beneath patients. Gel products are cost effective, latex free, fast & easy cleaning.





2		and the second		
-		SK-	29	
	571	1	100	

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

5

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549 Sung Kwang Meditech Co.,Ltd

14

1


5

6

8

www.skmeditech.com

Sung Kwang Meditech $C \in$





16

Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

Sung Kwang Meditech $C \in$



Gel body pad provides complete

SK-23-2	84(L) * 50(W) * 1.3(H)cm	
SK-23	160(L) * 50(W) * 1.3(H)cm	
SK-23-1	183(L) * 50(W) * 1.3(H)cm	

This set includes head, torso and foot section pads for use with segment tables

SK-23-1-1 Gel body pad 3-piece set(a+b+c)		
SK-23-1-1a (Head)	50(L) * 23(W) * 1.3(H)cm	
SK-23-1-1b(torso)	100(L) * 50(W) * 1.3(H)cm	
SK-23-1-1c(foot)	50(L) * 50(W) * 1.3(H)cm	

It provides gentle pressure for neonatal patients in the NIC units. It is designed for the isolette beds.

SK-23-5	60(L) * 33(W) * 1.3(H)cm	
SK-23-5-1	60(L) * 40(W) * 1.3(H)cm	


Gel chair pad is used underneath prominence area in a variety of applications.

SK-24	42(L) * 42(W) * 1.3(H)cm
SK-24-2	38(L) * 38(W) * 1.3(H)cm
SK-24-2-1	40(L) * 40(W) * 1(H)cm
SK-24-2-2	40(L) * 45(W) * 1(H)cm
SK-24-3	32(L) * 32(W) * 1.3(H)cm

Gel embo chair pad is used for iliac crest, under heels and elbows. Is used as kidney rest, axillary roll while patient in lateral position.

Gel embo body pad provides complete pressure reduction and shears cross entiret ablesur face.

Gel embo cha	ir pad	
SK-24-1	40(L) * 45(W) * 1.6(H)cm	
Gel embo boo	ly pad	
SK-23-3	160(L) * 45(W) * 1.6(H)cm	
Gel embo hal	f body pad	
SK-23-4	80(L) * 45(W) * 1.6(H)cm	



Tel: 82-31-855-7548 Fax: 82-31-855-7549

Operating Room



18


Sung Kwang Meditech $C \in$



Operating Room

Sung Kwang Meditech Co.,Ltd

SILICON




Sung Kwang Meditech $C \in$

www.skmeditech.com



DIRECTIONAL-STROKING PERCUSSOR



WRIST BAND **Elbow Band**

www.skmeditech.com

Sung Kwang Meditech Co.,Ltd

The fully equipped patent transport board is ideal for transport when patients need to be transferred from MRI, operating table, bed, etc.

SIZE:	
SK-MTB80	80 x 35cm
SK-MTB140	140 x 45cm
SK-MTB183	183 x 55cm
SK-MTSR860	86 x 50cm (Unfolded)
SK-MTSR1260	126 x 50cm
SK-MTSR1800	180 x 50cm

SK-MSG-NEOCUSSOR is crafted for gentle, effective chest therapy in infants, providing safe and specialized techniques like vibration, percussion, and postural drainage. These devices aid in clearing lung secretions, improving respiratory function, and reducing infection risk, supporting overall lung health in neonates.

CUSHION PAD

Multi-Purpose position pillow for hospital operating room. It is easy to use for short-term surgery like

Hemorrhoid surgery

Inner : soft cotton cushion



40.5(L) * 17(H) * 15.5(W)cm Color : White SK-13 SK-14 200(L) * 11(W)cm SK-14-1 70(L) * 7(W)cm SK-CU-CP-L 45(L) * 70(W) cm SK-CU-CP-S 40(L) * 60(W) cm SK-14-6 18(L) * 18.5(W) cm SK-14-7 35(L) * 12(W) cm Custom made available

Contact to E-mail: skmeditech@skmeditech.com


Tel: 82-31-855-7548 Fax: 82-31-855-7549

Rehabilitation

GEL WHEELCHAIR CUSHION PAD WITH COVER

It was made with a high-quality polyurethan gel and PU foam. The inner gel material is highly resistant to bacteria and bodily fluids, and can be easily cleaned with warm soapy water or a disinfectant. Also, it helps dissipate heat buildup.

Gel Wheelchair cushion pad offers excellent pressure reduction & distribution, great comfort, and positioning stability. It can be applied for individuals with a history of skin breakdown or those who suffer from pressure sores. Feel free to enjoy our products below.



Sung Kwang Meditech $C \in$

GEL WHEELCHAIR CUSHION 2" GEL WHEELCHAIR CUSHION 3"

SK-1616-2	16(L) * 16(W) * 2(H) Inch
SK-1816-2	18(L) * 16(W) * 2(H) Inch
SK-1818-2	18(L) * 18(W) * 2(H) Inch
SK-1616-3	16(L) * 16(W) * 3(H) Inch
SK-1816-3	18(L) * 16(W) * 3(H) Inch
SK-1818-3	18(L) * 18(W) * 3(H) Inch

SK-C1616-3.5	16(L) * 16(W) * 3.5(H) Inch
SK-C1816-3.5	18(L) * 16(W) * 3.5(H) Inch
SK-C1818-3.5	18(L) * 18(W) * 3.5(H) Inch
SK-MC1616-4	16(L) * 16(W) * 4(H) Inch
SK-MC1816-4	18(L) * 16(W) * 4(H) Inch
SK-MC1818-4	18(L) * 18(W) * 4(H) Inch

SK-E1616-2	16(L) * 16(W) * 2(H) Inch
SK-E1816-2	18(L) *16(W) *2(H) Inch
SK-E1818-2	18(L) * 18(W) * 2(H) Inch
SK-E1616-3	16(L) * 16(W) * 3(H) Inch
SK-E1816-3	18(L) * 16(W) * 3(H) Inch
SK-E1818-3	18(L) * 18(W) * 3(H) Inch

	and the second se	
	the second s	
	COMPANY OF THE PART CAN'T	
	and the second	
	and the same start of the same	
	and the second	
	and the second	
	and the second	
	A CONTRACTOR OF	
	and the second state of th	
	CONTRACTOR OF A	
	the second s	
	and the second	
	and the second	
	and the second se	
	a second s	
	and the second se	
and the second se	and the second	
	the second se	
100		
	100 million (1997)	

SK-EC1616-3.5	16(L) * 16(W) * 3.5(H) Inch
SK-EC1816-3.5	18(L) * 16(W) * 3.5(H) Inch
SK-EC1818-3.5	18(L) * 18(W) * 3.5(H) Inch
SK-EMC1616-4	16(L) * 16(W) * 4(H) Inch
SK-EMC1816-4	18(L) * 16(W) * 4(H) Inch
SK-FMC1818-4	18(L) * 18(W) * 4(H) Inch

GEL WHEELCHAIR BACK CUSHIONS WITH COVER

Polyurethane gel and PU foam helps comfort and durability for long-term.

The embossed gel feature aids comfort for long-lasting back and lumbar support in long-term seating environments.



RH-6-1

SNOW PAC



GEL HORSE S EQUESTRIAN PAD


RECTANGULAR ROUND



SILVER PRODUTS

GEL EMBO CUSHION PILLOW / MATTRESS WITH COVER (Memory foam + embo gel)

Characteristics: C-spine supporting effect is developed by medical basis. As humanengineering structure, it softly supports head and provides comfortable.



Sung Kwang Meditech Co.,Ltd	www.skmeditech.com	Contact to E-mail: skmeditech@skmeditech.com	Sung Kwang Meditech Co.,Ltd	www.skmeditech.co
24		Tel: 82-31-855-7548 Fax: 82-31-855-7549		

RH-4(Dual)+cover 39(L) * 10(W) * 1(H)cm

RH-6 (close type) 45(L) * 45(W) cm RH-6-1(open type) 45(L) * 46(W) cm



RH-4

RH-6

Sung Kwang Meditech $C \in$

It is used for lugeric patient and unpositioned handicapped person. Make position easy and fast. Then, suction equipment is used.

RH B type	83(L) * 108(W)cm
HB-102a	33cm



-EQ-a Rectangula	ar shape with cover
36	(L) * 50(W) * 1(H) cm
-EQ-a cover Rect	angular shape cover only
36	(L) * 50(W) * 1(H) cm
-EQ-b Round sha	pe with cover
40	(L) * 55(W) * 1(H) cm
-EQ-b cover Rour	nd shape cover only
40	(L) * 55(W) * 1(H) cm

Outer : 100% cotton velvet with silvercholoid coating Inner : cotton 75% polyester 25% Materials : slowly recovery polyurethane foam 100% + embo gel pad **Color** : silver, ivory Kind : therapeutic, antibacterial, deodorizational.

RH-1	50(L) * 33(W)cm
RH-2 (twin)	RH-2-1 (full)
RH-2-3 (queen)	RH-2-4 (king)

Cover Only

,	
RH-2-cover	RH-2-1-cover
RH-2-3-cover	RH-2-4-cover

com

Lead Aprons

Several Types of Aprons For Radiation Protection



SK-15 (Buckle type)



SK-15-1 (Coat type)



SK-15-4 (X-band type)



SK-15-5 (Buckle X-band type)

Exposure to radiation is fatal to the medical professional and the human body.

Here are all the radiation protection products you need for protection from radiation exposures.

We manufacture a comprehensive range of high-quality, lightweight, and cost-effective products for medical customers around the world.

- Inner : Natural rubber-based Lead flexible wearing comfort and best.
- Outer : nylon fabric, etc satisfies the workmanship hygienic requirements and the demands of medical use.

Pb Equivalent:

0.25 / 0.35 / 0.50 mmPb Available

±2~4cm, unit∶cm L : Length, W : width

	S		М		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	93	58	98	60	103	63	108	66
Coat type (Front/back)	96/86	53	99/89	55	105/95	57	106/96	60
X-band type	94	55	96	57	100	60	104	62
Vest/Skirt type	59/50	49/42	60/52	51/43	61/53	53/46	63/54	56/49
panorama type	82	47	85	50	88	53	91	55

FEMALE

MALE

	S		М		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	84	50	89	52	95	56	99	58
Coat type (Front/back)	85/75	47	91/81	49	95/85	52	100/90	54
X-band type	85	51	90	53	96	56	100	60
Vest/Skirt type	55/46	46/38	56/48	48/40	57/49	51/43	59/50	54/46
panorama type	76	43	81	47	84	50	87	53

DROP OFF



Remove your apron quickly and easily with shoulder valcro.


Available to SK-15, SK-15-1, SK-15-4, SK-15-5 Several Types of Aprons For Radiation Protection






















Sung Kwang Meditech Co., Ltd

www.skmeditech.com

SK-15-6 (X-buckle type)

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549

26-

Sung Kwang Meditech $C \in$

Non Lead Aprons

Several Kinds of X-Ray Apron For Surgical Patient



SKN-15 (Buckle type)



SKN-15-1(Coat type)



SKN-15-4 (X-band type)



SKN-15-5 (Buckle X-band type)

Exposure to radiation is fatal to the medical professional and the human body.

Here are all the radiation protection Lead-Free products you need for protection from radiation exposures.

We manufacture a comprehensive range of high-quality, lightweight, and cost-effective products for medical customers around the world.

- **Inner** : Natural rubber-based Lead flexible wearing comfort and best.
- **Outer** : Water Resistant, Abrasion Resistant, Cool wear

Pb Equivalent:

0.25 / 0.35 / 0.50 mmPb Available

2.4

	MALE						cm, unit : cm th, W : width	
	S		M		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	93	58	98	60	103	63	108	66
Coat type (Front/back)	96/86	53	99/89	55	105/95	57	106/96	60
X-band type	94	55	96	57	100	60	104	62
Vest/Skirt type	59/50	49/42	60/52	51/43	61/53	53/46	63/54	56/49
panorama type	82	47	85	50	88	53	91	55

FEMALE

	S		М		L		XL	
	L	W	L	W	L	W	L	W
Buckle type	84	50	89	52	95	56	99	58
Coat type (Front/back)	85/75	47	91/81	49	95/85	52	100/90	54
X-band type	85	51	90	53	96	56	100	60
Vest/Skirt type	55/46	46/38	56/48	48/40	57/49	51/43	59/50	54/46
panorama type	76	43	81	47	84	50	87	53

DROP OFF



Remove your apron quickly and easily with shoulder valcro.

Available to SKN-15, SKN-15-1, SKN-15-4, SKN-15-5



SKN-15-2 (Vest & Skirt type)

Several Kinds of X-Ray Apron For Surgical Patient



SKN-15-3 (Half type) SKN-15-3-1 (Panorama type) SKN-15-3-2









Sung Kwang Meditech Co.,Ltd

www.skmeditech.com

SKN-15-6 (X-buckle type)

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549 Sung Kwang Meditech Co.,Ltd

Sung Kwang Meditech $C \in$









Thyroid shield is an essential to protect the thyroid gland from radiation exposure. (SKN-16, SKN-16-a, SKN-16-b)

Find a wide variety of radiation protection products and colors.

Pb Equivalent : 0.25 / 0.35 / 0.50 mmPb Available

Color type



Dia Navy

Pearl Gre



Pearl Green





Black Dot



Dia Orange







Pearl Gold

Melange Orange

Honeycomb Grey



Dia Chocolate ****

		5.	6	
			Ġ.	0.
2	9		ā.	
2				

Pearl Purple



Melange Navy



Our Lead-Free cap and gonad shield are lightweight and comfortable.

An adjustable velcro in the back is attached for a customizable fit.

- -Two types for cap : SKN-16-4 / 16-4-a
- -Two types for gonad : SKN-16-5 / 16-5-a
- -Order Types : Male, female
- -Protection : 0.25mm front
- -Pb Equivalent: 0.25 / 0.35 / 0.50 mmPb Available

Custom-made available.

Lead Free cap	SKN-16-4 (Basic Model)	
Lead Free cap	SKN-16-4-a (Cool wear mositure-wicking fabr	ric)
Lead Free gona	ad shield SKN-16-5	
Lead Free gona	ad shield SKN-16-5-a	





Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549 _29

www.skmeditech.com

Radiation Protection Accessories

APRON HOLDER APRON WALL CONSOLE SK-16-6-1 SK-16-6-2 SK-16-6 SK-16-6-a SK-16-6-b For Vest & Skirt Type FOR PREGNANT TABLE LEAD / WOMAN BARRIER NON-LEAD CURTAIN drawing. Order made Order made SK-16-7-3 SK-16-7-4 (양면) SK-16-7-4-a (단면)

LEAD BARRIER



Apron holder will accommodate up to 2-10 aprons. Unit is ideal if wall space is limited. It is designed to prevent declining on side in moving. Solid stainless steel construction prevents rust in long-term use. High durability under heavy use. Can be assembled after delivering.

up to 10
up to 5
up to 10
up to 2
up to 5

Lead curtain provides separated radiation shielded space by simply

- protection : 0.5 mmpb
- Protection for pregnant woman barrier : 0,5 mmpb

SK-16-6	71(L) * 63(W) * 140(H) cm
SK-16-6-a	50(L) * 63(W) * 140(H) cm
SK-16-6-1	20(L) * 50(W) * 25(H) cm
SK-16-6-2	25(L) * 50(W) * 25(H) cm
SK-16-7-1	57(W) * 60(L)/83(H) to 146(H) cm
SK-16-7-3	200(W) * 150(H) cm
SK-16-7-4	100(L) * 65(W) * 70(H) cm
SK-16-7-4-a	100(L) * 65(W) * 70(H) cm

Lead-glass viewing panel provides exceptional visibility.

- Protection : 1,76 mmpb
- Order made

cm
cm

Contact to E-mail: skmeditech@skmeditech.com

Tel: 82-31-855-7548 Fax: 82-31-855-7549

DISPOSABLE COLLAR & CAP RADIATION PROTECTION

Thyroid collars and caps are available as disposable products, made from soft, non-woven fabric.

These items offer a hygienic, single-use option with absorbent, non-irritating liners that protect against dirt and perspiration.



X-RAY GLOVES



These gloves are designed for use when exposed directly to X-rays.

They are suitable for use in veterinary clinics, nuclear facilities, and other similar environments.

LEAD RULER

Ruler used for medical imaging to track the location of the affected area and measure its length when interpreting medical photos and videos.

Size : 45, 60, 115, 150, 180 cm Available

LEAD MARKER Lead markers are used in medical imaging, such as X-rays and CT scans, to indicate the exact location of the affected area and assist in accurate diagnosis and treatment planning. Available : Blue, Red, Reverse, Numbering

RADIATION SHIELDING TAPE

Made from non-hazardous, lead-free materials, Radiation Shielding Tape is designed to Key Features effectively block or reduce radiation exposure in hazardous environments. It is ideal for use in medical, industrial, and laboratory settings, offering reliable protection for both personnel and sensitive equipment.



Sung Kwang Meditech Co., Ltd

www.skmeditech.com

Sung Kwang Meditech Co.,Ltd

Sung Kwang Meditech $C \in$

Advantages

- 1. Disposable—no need for cleaning
- 2. Convenient and cost-effective
- 3. Lead-free and eco-friendly
- 4. Lightweight
- 5. Available in multiple sizes

Pb Equivalent Available

- 0.25 mmPb
- 0.35 mmPb
- 0.50 mmPb

Fabric Color : Blue

Size : One size fits all

Packaging Options

- Box of 10
- Box of 20
- Box of 50

RADIATION SHIELDING X-RAY CURTAINS FOR SCANNER



Radiation Shielding Curtains for security, food scanners. The protection range includes lead equivalent values of 0.125, 0.25, 0.375, 0.5mm depending on the level of protection required.

Customer made.



- Lead-Free, Non-Hazardous Material 3 Ensures safety while maintaining high shielding efficiency.
- Effective Radiation Blocking : Provides strong protection against various types of radiation exposure.
- Broad Application : Suitable for use in healthcare, industrial, and research settings where radiation is present.

Scatter Radiation Protection Shields for Use in Interventional Procedures

ANGIOGRAPHY SHIELD WITH OPENING



Provides optimal protection for femoral artery access procedures in IR, OR, Cath, and EP labs, without obstructing the physician's technique.

- Size : 370mm x 420mm / Customizable
- Lead Equivalency (mm Pb) 0.25mm / 0.375mm
- Attenuation : Over 90% at 90 kVp
- 100% lead-free, easy to dispose
- High absorbency covering material
- Application : Angiography, Coronary Catheterization

BILIARY SPLIT SHIELD



This versatile shield delivers exceptional protection for both AP and oblique viewing angles.

- Size : 370mm x 420mm / Customizable
- Lead Equivalency (mm Pb) 0.25mm / 0.375mm
- Attenuation : Over 90% at 90 kVp
- 100% lead-free, easy to dispose
- High absorbency covering material
- Application : Biopsy Procedures, Transjugular Intrahepatic Portosystemic Shunt (TIPS)

PERIPHERAL SHIELD



Designed to offer comprehensive protection during peripheral procedures. The extended length ensures full mobility and complete coverage when imaging the limb.

- Size : 350mm x 900mm / Customizable
- Lead Equivalency (mm Pb) 0.25mm / 0.375mm
- Attenuation : Over 90% at 90 kVp
- 100% lead-free, easy to dispose
- High absorbency covering material
- Application : Peripheral Procedures, Access through Femoral Artery






•


EP SUBCLAVIAN SHIELD



BIOPSY SHIELD



DIALYSIS SHIELD



Sung Kwang Meditech Co., Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com Tel: 82-31-855-7548 Fax: 82-31-855-7549 Sung Kwang Meditech Co.,Ltd

32-

www.skmeditech.com

MULTIPURPOSE SHIELD



An ideal alternative or supplementary shield for use in IR, Cath, and EP labs, tailored for procedure-specific needs.

- Size : 350mm x 450mm / Customizable
- Lead Equivalency (mm Pb) : 0.25mm / 0.375mm
- Attenuation : Over 90% at 90 kVp
- 100% lead-free, easy to dispose
- High absorbency covering material
- Application : Left and Right Radial, Shunt, Biliary Drainage, Spinal, Fibroid Embolization

Offers superior protection during device implant procedures like Bi-V pacing. Specially designed opening ensures easy access and enhanced protection.

- Size : 280mm x 410mm / Customizable
- Lead Equivalency (mm Pb) : 0.25mm / 0.375mm
- Attenuation : Over 90% at 90 kVp
- 100% lead-free, easy to dispose
- High absorbency covering material
- Application : Bi-Ventricular Pacemaker Implants, ICS Procedures, Bi-Ventricular Pacing



Delivers excellent protection for all biopsy procedures. The rectangular opening allows clear access and visibility while reducing radiation exposure.

- Size : 350mm x 450mm / Customizable
- Lead Equivalency (mm Pb) : 0.25mm / 0.375mm
- Attenuation : Over 90% at 90 kVp
- 100% lead-free, easy to dispose
- High absorbency covering material
- Application : Needle Biopsy or Puncture Procedures

Adjustable shield that provides superior protection, with added flexibility for oblique fluoroscopy views.

Lead Acrylic Radiation Protection Products

Since human's eyes are more sensitive to radiation than other parts of human body, using eyewears or masks is strongly recommended for users related to the following vocations: Fluoroscopic Procedures, Cardiac Cath, EP Labs, Urology, Interventional Radiology, Orthopedic Surgery.

Care Instructions:

- Wipe clean with a soft cloth.
- Avoid using abrasive materials, as they may scratch the front lens.
- Not suitable for sterilization with high temperatures or high pressure.

This mask offers enhanced protection for the eyes and face during radiological procedures, minimizing radiation exposure. With a multi-adjustable headpiece, the mask can be angled to shield from radiation in various positions and is easily adjustable for a comfortable, personalized fit.

SK-AC-FFM (Full Face Mask)

- Weight: 400g
- Protection: 0.15mmPb
- For Full Face (Downward Focus)



SK-AC-HFM (Half Face Mask)

- Weight: 400g
- Protection: 0.15mmPb
- For Half Face (Half face and half sides Focus)




Made in Korea

SK-AC-BASIC

- Weight: 60g
- Protection: 0.08mmPb



SK-AC-FLIP





SK-AC-BASIC, SK-AC-FLIP:

X-ray protective lightweight Acrylic goggles offer comfort and provide full protection from the front and side, specifically designed to decrease scattered radiation.

Mainly for those without glasses, but possible for those with glasses (Depending on Glasses sizes).

Tel: 82-31-855-7548 Fax: 82-31-855-7549





Made in Korea

Sung Kwang Meditech Co., Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com

SK-400PSM (PANORAMIC FACE SHIELD)



◆ Weight: 350 grams ◆ Protection : .10mm Lead Acrylic

SK-400PSM

SK-450FPM (FULL FACE SHIELD)



SK-450FPM

SK-480 (With G

- Weight: 60g
- protection: 0,07r



SK-380 (Without Glasses)

• protection: 0.07mmPb

• Weight: 60g





Import Sung Kwang Meditech C€

The 400PSM Panorama Shield consists of a Face Formed Radiation Protective Lead Impregnated Acrylic Shield that is suspended from a dual adjustable headpiece.

The acrylic shield can be positioned at any angle to the face so as to minimize scatter radiation from passing under the shield by side adjustable dials. The headgear is adjustable in two positions.

An adjustable headband that saddles the crown of the head while connecting at a 90 angle to an adjustable headband that surrounds the head circumference from brow to back. The circumference portion of the headband has an adjustable dial in the rear to increase or decrease the hat size for a comfortable, snug fit.

Wipe clean with a soft cloth. Abrasive materials may scratch front lens.

Cannot be sterilized with high temperatures or high pressure.

◆ Material : Acrylic/Nylon

The 450FPM Panorama Mask consists of a Face Formed Radiation Protective Lead Impregnated Acrylic Shield that is suspended from an adjustable headpiece.

The stationary acrylic shield extends from a headpiece brow connector for a length of 8 inches from top to bottom in front of the entire face. A support strap saddles the crown of the head while connecting at a 90 angle to an adjustable headband that surrounds the head circumference from brow to back.

The circumference portion of the headband has an adjustable strap in the rear to increase or decrease the hat size for a comfortable, snug fit. W ipe clean with a soft cloth.

Abrasive materials may scratch front lens. Cannot be sterilized with high temperatures or high pressure.

◆ Weight: 450 grams ◆ Protection: 10mm Lead Acrylic Material: Acrylic/Nylon

ilasses)	
mmPb	Curved lead acrylic lenses effectively block scattered radiation from multiple angles.
	For optimal protection, position the lens center at eye level and wear it close to the face.
Tour	7
	Made in Japan
om	Contact to E-mail: skmeditech@skmeditech.com
	Tel: 82-31-855-7548 Fax: 82-31-855-7549 35

Import Lead Eyewear (With Glasses)

SK-RE-FGUARD (Beta/Alpha)



The Fitguard fitover is a versatile, high quality nylon frame. Featuring Flexfit technology, the Fitguard lead glasses' temples can be benched (adjusted) easily for a secure, comfortable fit. Blacked out side shields provide a streamlined look as well as 0.50mm lead equivalent protection

The Beta model is ideal for fitting over standard prescription eyewear with square or rounded lenses.

The Alpha model is ideal for fitting over wide prescription eyewear while still providing a secure and comfortable fit. With the most balanced weight distribution, the Fitguard Alpha is our most popular fitover model.

- ◆ Weight: 108g(B) / 91g(A)
- ◆ Front Lens: 0.75mm Pb Glass
- ◆ Side Shields: 0,50mm Pb Strips
- ◆ Frame Material: Nylon
- ◆ Frame Szie: 64-18-150(B) / 69-17-145(A)
- ◆ Color: Black
- ◆ CE: Yes

SK-90 (90 FITOVER)



Designed to fit over prescription eyewear, model 90 fitovers are an excellent choice for those requiring protection against x-ray and liquid splashes.

The large, clear side shields protect the user from lateral exposure while eliminating the "tunnel" effects of look-alikes.

- ♦ Weight: 112g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb Standard
- ◆ Frame Material: Durable Nylon
- ◆ Frame Size: 66-14-120
- ♦ Color: Matte Black
- ◆ CE: Yes

SK-70S (70S ASTRO FLEX) (Without Glasses)





A panoramic goggle constructed of impact resistant propionate.

The frame front has molded side shielding and is fitted to the brow for increased splash protection.

The temple length and frame angle are adjustable for a uniform custom fit

Model 70S comes with Lead side shields for necessary protection from oblique angle exposure.

- ◆ Weight: 88g
- ◆ Front Lens: 0.75mm Pb Glass
- Side Shields: 0.50mm Pb Strips
- ◆ Frame Material: Propionate
- ◆ Frame Szie: 54-19-Adj
- ◆ Color: Clear/Blue
- ◆ CE: Yes

Lead Eyewear (Without Glasses)

SK-RE-9941A (9941 ULTRALITE)



SK-N-BR (NIKE BRAZEN)





SK-RE-RAZER (RAZER)





	•	1
	٠	(
	٠	ŀ
	٠	F
	٠	(
	٠	(

Sung Kwang Meditech Co.,Ltd	www.skmeditech.com	Contact to E-mail: skmeditech@skmeditech.com	Sung Kwang Meditech Co.,Ltd	www.skmeditech.co
6		Tel: 82-31-855-7548 Fax: 82-31-855-7549		

Import Sung Kwang Meditech C€

- The 9941A Ultralites come standard with a rubberized nose bridge. They are designed to offer comfort and provide protection for the entire eye area without the need of additional lead side shields. The 9941A is a popular frame and ideal for those who have a high
- power prescription.
- ♦ Weight: 64g
- ◆ Front Lens: 0.75mm Pb Glass
- ◆ Side Shields: 0.50mmPb Strips Optional
- ◆ Frame Material: Durable Nylon
- ◆ Frame Szie: 55-22-130
- ◆ Color: Silver
- ◆ CE: Yes

The Nike Brazen features an advanced ventilated nose bridge for grip, comfort, and reduced fogging

- Its lightweight nylon design combine streamlined style and enhanced comfort
- The premium Nike Brazen frame offers versatility with its superior coverage and protection.
- ◆ Weight: 75g
- ◆ Front Lens: 0.75mmPb Glass
- ◆ Side Shields: 0.5mmPb Optional
- ◆ Frame Material: Durable Nylon
- ◆ Frame Size: 60-16-140
- ◆ Color: Matte Black
- ◆ CE: Yes

- Features: Ventilated Nose Bridge, Rubberized Temples for Grip and Comfort, Modern Matte Color Finishes, High Rx Options
- ♦ Weight: 62.5g
- ◆ Front Lens: 0.75mmPb Glass
- Side Shields: 0.5mmPb attached
- Frame Material: Polycarbonate
- Frame Size: 58-17-135
- Color: Matte Green, Matte Gray
- CE: Yes

com	Contact to E-mail: skmeditech@skmeditech		
	Tel: 82-31-855-7548	Fax: 82-31-855-7549	27

The Razer's blends a sleek design into a lightweight and durable frame. Its premium features such as a ventilated nose bridge, rubberized temples, and wide-angle lenses maximize comfort and protection from x-ray.

The modern style and use of light-weight materials make the Razer an excellent choice for longprocedures.

Import **Non Lead Aprons**

Our lighter lead-free material (0.25 / 0.35 / 0.5) Lead Free ◆ Environmentally Safe

ESSENTIAL VEST SKIRT APRON



Be prepared for prolonged procedures with the Vest Skirt apron.

The two piece design shifts 70% of the weight directly to the hips, reducing back fatigue.

The Essential apron configuration features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection.

The skirt comes standard with a premium non-slip fabric in the back. Essential aprons are only available in standard sizes.

◆ Protection : 0.25 / 0.35 / 0.5 mmPb

ESSENTIAL VEST SKIRT APRON

SK-AP-VSQ

ESSENTIAL WRAPAROUND BACK RELIEF APRON



BACK RELIEF APRON



Offering full protection, the Essential wraparound lead apron's backrelieving design transfers the weight from your shoulders to your hips with the aid of a 6" lumbar support belt.

An Essential apron configuration features extended front panels to achieve a 3/4 overlap (where protection is doubled) which allows for a lighter weight apron without compromising protection.

The Essential Wraparound is available in standard sizing only.

◆ Protection : 0.25 / 0.35 / 0.5 mmPb

ESSENTIAL WRAPAROUND BACK RELIEF APRON

SK-AP-WAQ

SK-AP-BR

The Back Relief apron features a padded belt system constructed with an elasticized panel that will not lose its shape.

When engaged, the back relief belt provides lumbar support while transferring the weight of the apron from your shoulders to your hips for increased comfort and back support.

A velcro backstrap is included and secures the apron around your shoulder area to prevent sliding.

Protection : 0.25 / 0.35 / 0.5 mmPb

BACK RELIEF APRON

Sung

FLEXBACK



area to prevent sliding.

FLEXBACK APRC

BUCKLE APRON

support.

BUCKLE APRON



TRI-TAB VEST SKIRT APRON



UROLOGY APRON



UROLOGY APRON

Kwang Meditech Co.,Ltd	www.skmeditech.com	Contact to E-mail: skmeditech@skmeditech.com	Sung Kwang Meditech Co.,Ltd	www.skmeditech
		Tel: 82-31-855-7548 Fax: 82-31-855-7549		

Import Sung Kwang Meditech CE

- Representing the gold standard of guality, the Flexback apron design provides enhanced back relief, comfort, and excellent weight distribution.
- Criss-cross back straps are specially designed with elastic to "flex" the apron's weight burden off your shoulders to your hips, thereby providing support to your lower back.
- A velcro backstrap is included and secures the apron around your shoulder
- ◆ Protection : 0.25 / 0.35 / 0.5 mmPb

N	SK-AP-FB

- The Buckle Apron comes standard with an easy-release buckle allowing for quick removal of the apron.
- The belt is adjustable to accommodate different sizes.
- Criss-cross back straps provide a secure, comfortable fit.
- ◆ Protection : 0.25 / 0.35 / 0.5 mmPb

- With an easy-on, easy-off multi closure design, the Tri-tab is a popular style for vest and skirt combinations
- Velcro closures are located on the left shoulder and around the hips.
- An elastic back allows for a flexible fit while providing comfort and back
- The Tri-tab Vest Skirt reduces excess overlapping lead and is up to 10% lighter than traditional vest models.
- Protection: 0.25 / 0.35 / 0.5 mmPb

|--|

SK-AP-TVS

- If you sit while performing procedures requiring x-ray, you may be unnecessarily exposing yourself to radiation.
- The Urology Apron features a three panel system at the bottom of the apron which enables the user to sit comfortably while maintaining maximum coverage from x-ray.
- ◆ Protection : 0.25 / 0.35 / 0.5 mmPb

SK-AP-UA

Import **Non Lead Aprons**

SURGICAL DROP APRON



TIE APRON



VELCRO ADJUSTABLE



MATERNITY APRON



Maintain your sterile field with the guick release Surgical Drop Apron.

Easy-release velcro closures allow the surgical drop to fall even when underneath your scrubs.

With padded shoulders and a criss-cross back strap design the surgical drop apron provides maximum comfort for all users.

Hanging straps are included to allow for easy storage on any Protech apron rack.

SURGICAL DROP APRON

SK-AP-SD

SK-AP-EA

The Tie Apron is an easy-to-use, affordable option for x-ray protection. Utilizing a tie closure, the Tie Apron's criss crossing straps ensure a secure and comfortable fit.

Do you need comfort and adjustability at an affordable price?

The Velcro Adjustable radiation protective apron has crisscrossing back straps that provide back support while securing the apron around your waist. Like the Flexback apron, wide back straps prevent the apron from slipping

around your shoulder and back area.

VELCRO ADJUSTABLE APRON

TIE APRON

SK-AP-VA

The maternity apron features easy ties and a guick release buckle to secure the back and front panels; allowing complete adjustability during pregnancy.

The maternity apron also includes an additional .50mm equivalent half apron for a total of 1.00mm Pb protection over the abdominal area and lower extremities.

Please note, because this apron contains a lot of fabric, we do not recommend it in Standard Lead as it would be much too heavy.

Maternity Apron

SK-AP-MA

Sung Kwang Meditech Co., Ltd

www.skmeditech.com

Contact to E-mail: skmeditech@skmeditech.com

Lead Gloves and Sleeves



SKN-16-1 (Size: 6~9)



SK-SL (Size: S/M/L)

Sung Kwang Meditech Co.,Ltd



XGUARD Radiation Reducing Sleeve, when used in conjunction with the XGUARD Radiation Reducing Procedure Glove creates an attenuative barrier from your fingertips to your biceps from scatter radiation.

- ATTENUA

XGUARD 60 Kv Sleeve 55%

SURGICAL DROP THYROID CO BREAST SHIELD **SHOULDER** MODEL

www.skmeditech.c

Tel: 82-31-855-7548 Fax: 82-31-855-7549

40

Our Soft-Touch finish creates a smooth interior reducing skin irritation. Anatomically correct curved fingers reduce hand fatigue.

Sizing range from 6 to 9 in half sizes for an unequaled comfortable fit. Attenuates up to 50% more direct beam radiation than other leading products. ThInner fingertips provide enhanced tactile sensitivity.

ATTENUATION

RR-2:0.012"/0.30mm					RR-3:0.024" / 0.60 mm			
Кур	60	90	120		Kvp	60	90	120
Cuff	56%	44%	37%		Cuff	76%	63%	55%
Palm	57%	45%	38%		Palm	75%	62%	54%
Finger	57%	56%	38%	ĺ	Finger	75%	62%	54%

ATTENUATION

MODEL 1 : 0.009" / 0.22 mm				MODEL 2 : 0.012" / 0.30 mm			
Kvp	60	90	120	Kvp	60	90	120
Cuff	51%	40%	34%	Cuff	62%	51%	43%
Palm	55%	44%	37%	Palm	63%	51%	43%
Finger	55%	44%	37%	Finger	61%	49%	42%

• Inner Nitrile coating - reducing proteins, antigens and irritation.


• Lightweight material reduces arm fatigue.

• Extreme flexibility does not restrict movement or dexterity.

Available in three sizes for a comfortable fit.

TION						
γp	80 Kvp	100 Kvp	120 Kvp	THICKNESS		
)	45%	35%	31%	0.012" - 0.30mm		


OLLARS		
1		
	BREAST SHIELD	APSTOLE
9	SURGICAL DROP SHOULDER	APDRSH
	THYROID COLLERS	APTCO / APTCS
com	Contact to E-mail: sk	meditech@skmeditech.com
	Tel: 82-31-855-7548	Fax: 82-31-855-7549 41



Import Sung Kwang Meditech ⊂ €


ier.				
				420 Denier
	Hot Pink #125	Marine Blue #135	Navy #140	Neon Green #145
	Royal #175	Tan #180	Teol #185	420 Denier Yellow #190
	Green #225	Navy #230	Purple #235	Red #240
ΥV	'INYLS*	Black #300	Silver #320	Burgundy Starlight #305
ver	Blue Flame #425	Brown Animal Spots #435	Butterfly #40	Koo Koo Kano - Backwoods #443
nks	Camo - Yellow #480	Chicklets Blue #490	Crayon Hearts #500	Divebar #502
i i	Galaxy - Royal #525	Green Neon	Hawaiian #545	Hawaii Postcard #550
	Mardi Gras #590	Tie-Dye #540 Wight of the second seco	#949 Painted Dragon #600	Paint Splash #805
le	Plaid - Tan	Plaid - Teal	Purple Neon	Reinbow
S Internet	#635 For Gazer #685	#640	Tie-Dye #645 Steel #695	#650 Surf Wex #705

Rose #170





Sung Kwang Meditech Co.,Ltd





T) 82-31-855-7548 M) 82-10-4590-7548 F) 82-31-855-7549 E) skmeditech@skmeditech.com Address) 1014-20, Bukhansan-ro, Jangheung-myeon, Yangju-si, Gyeonggi-do, Korea Sung Kwang Meditech co.,Itd

www.skmeditech.com